
Uploading Images Using A Provider’s Webpage Uploader

Patient · https://mymedicalimages.com/help/uploading-images-using-a-providers-webpage-uploader/

Uploading Images Using A Provider’s Webpage Uploader


Step 1: Navigate to the provider’s mymedicalimages webpage uploader. This can be done by entering the address for the webpage. The address for the webpage will be https://app.mymedicalimages.com/fupi/(doctor’s extension).

Alternatively, your provider may have a "Upload on mymedicalimages" button, similar to one of the buttons displayed below. You can click on these to be taken directly to the web uploader.




Step 2: If you do not have an account, fill in your information in the text fields given, along with a brief description of the images you are uploading. Then click on the "I accept the mymedicalimages

Terms of Service" checkbox, and click on the blue "Upload" button.



If you already have an account with mymedicalimages, click on the blue "Click HERE" button. Enter in your account information, along with a brief description of the images you are uploading, then click on the blue "Upload" button.



Step 3: Select the file source.


You will have the option to pick between a CD, a .zip file, or "other". If the images you are attempting to upload are not on a CD or contained in a .zip file, you will need to select the "other" option.

The following steps you need to take will change depending on what source type you choose.

If your files are on a CD, select the "CD" option and go to "For CD".


If your files are on a .zip file, select the ".zip" option and go to "For .zip File".

If your files are on neither a CD or a .zip file, select the "other" option and go to "For Other".




For CD:

1. Insert the disk and click on the next button.




2. Click on the "Click Here" button to open the file explorer window.




3. Select your computer’s CD Drive.

4. Select the folder in your CD containing your DICOM files.



5. Click on the "upload" button in the bottom right-hand corner of the window.

6. In the pop-up window labeled "Upload files to this site?" click on the blue "Upload" button.




7. Go to Step 4.

For .zip:


1. Click on the "Click here" button to open the file explorer window.

2. Navigate to the .zip file you wish to upload and select it.

3.Click the "Open" button in the bottom right-hand corner of the window.




4. Go to Step 4.




For Other:

1. Click on the "Click here" button to open the file explorer window.




2. Navigate to the file you wish to upload and select it.

3. Click the "Upload" button in the bottom right-hand corner of the window.




4. In the pop-up window labeled "Upload files to this site?", click the blue "Upload" button.




5. Go to Step 4.

Step 4: You will be taken to the upload page. This page will display the amount of studies that are about to be uploaded inside of the progress bar in the center of the screen. Click on the green "Start" button to begin uploading your images.



Once the upload begins, the progress bar will begin to fill as the files are uploaded to mymedicalimages.com.

Depending on the amount of files being uploaded, this may take a while to complete. Please make sure that you leave the page open and remain connected to the internet for the entire duration of the upload, or you may encounter issues, such as corrupted files.



Once the upload is complete, you will be redirected to the thank you page. You can click on the home button to return to your mymedicalimages account, where a copy of your images can be found on the dashboard.

If you did not have an account before the upload, a new, basic account will be created for you with your images copied into it.

